

OneCheck
COVID-19 Ag

CROMATOGRAPHIC IMMUNOASSAY RAPID TEST

REF TR186

IVD



20 pcs

INTENDED USE

The COVID-19 Ag Rapid Test is a lateral flow chromatographic immunoassay for the qualitative detection of Sars-CoV-2 Antigen in human nasopharyngeal swab specimen as an aid in the rapid diagnosis in individuals with symptoms of SARS-CoV-2 infection.

Antigen is generally detectable in upper respiratory specimens during the acute phase of infection, but the titer may be too low to detect in the early stages or in asymptomatic individuals.

In case of positive result it is recommended to confirm it with an alternative method. Clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Negative results should be considered in the context of a patient's recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19.

This test is for professional and in vitro diagnostic use only and it is intended for use by trained clinical laboratory personnel.

SUMMARY

COVID-19 is an acute respiratory infectious disease caused by the novel coronavirus SARS-CoV-2. The incubation period of the disease is 1 to 14 days, mostly 3-7 days.

Currently, the patients infected by the novel coronavirus, symptomatic and asymptomatic, are the main source of infection. An infected person can be contagious even during the incubation period.

The main manifestations include fever, fatigue and dry cough, breathing difficulties. In more severe cases, infection can cause pneumonia, severe acute respiratory syndrome, kidney failure and even death.

IgM antibodies can be detected a few days after the incubation period and the appearance of the first symptoms. They remain present for a short time and can be an indicator of acute infection. IgG antibodies appear after several days and remain present for a long time.

The Antigen is present before the formation of antibodies and its detection is the most suitable test for an early diagnosis of the infection.

PRINCIPLE

The COVID-19 Antigen Rapid Test is a qualitative membrane-based immunoassay for the detection of SARS-CoV-2 Antigens in human nasopharyngeal swab specimen. During testing, the specimen reacts with SARS-CoV-2 antibody-coated particles in the test. The mixture then

migrates upward on the membrane by capillary action and reacts with the SARS-CoV-2 antibody in test line region. If the specimen contains SARS-CoV-2 Antigens, a colored line will appear in test line region as a result of this. If the specimen does not contain antigens to SARS-CoV-2, no colored line will appear in the test line region, indicating a negative result. To serve as a procedural control, a colored line will always appear in the control line region, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

MATERIALS

MATERIALS PROVIDED

- 20 test cassettes
- 20 sterile swabs
- 20 specimen extraction tubes
- 2 bottles with extraction buffer
- 1 workstation
- 1 user manual

MATERIALS REQUIRED BUT NOT PROVIDED

- Timer

PRECAUTIONS

- For professional in vitro diagnostic use only.
- Do not use after expiration date.
- Do not use test if pouch is damaged.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Handle all specimens as if they contain infectious agents.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- The used test should be discarded according to local regulations.


STORAGE AND STABILITY

- Store as packaged in the sealed pouch at room temperature (4-30°C). Do not freeze.
- The test is stable through the expiration date if it remain in the sealed pouch until use and properly stored.

SPECIMEN COLLECTION AND PREPARATION

SPECIMEN COLLECTION

1. Insert a sterile swab into the nostril of the patient, reaching the surface of the posterior nasopharynx.
2. Swab over the surface of the posterior nasopharynx.
3. Withdraw the sterile swab from the nasal cavity.




SPECIMEN PREPARATION

1. Add 10 drops of buffer in the specimen collection tube.
2. Insert the swab specimen into the specimen collection tube. Press against the inner wall of the tube and stir the swab for approximately

10 seconds while pressing the swab head against the inner wall of the tube to release the antigens in the collection tube.

3. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
4. Fit the tube tip on top of the extraction tube.




SPECIMEN STORAGE

Specimens should be tested as soon as possible after collection.

- If swabs are not been processed immediately, it is highly recommended the swab sample is placed into a dry, sterile, and tightly sealed plastic tube for storage. The swab specimen in dry and sterile condition is stable for up to 8 hours at room temperature and 24 hours at 2-8°C.
- The storage of the specimen after extraction is stable for 2 hours at room temperature or 24 hours at 2-8°C.

PROCEDURE

1. Allow the test, specimen, buffer and/or controls to reach room temperature (15-30°C) prior to testing.
2. Remove the test cassette from the foil pouch and place it on a clean and level surface.
3. Unscrew the dropper cap and transfer **3 drops of solution** (75 µl) to the specimen well avoiding the formation of air bubbles.
4. Wait **15 minutes** and read the result.
5. Do not interpret the result after 20 minutes.



INTERPRETATION OF RESULTS

One colored line appears in the control line region if the test was performed correctly.

NEGATIVE: Only one colored line appears in the control line region "C" and **no line appears in the test region "T"**. Antigen in the sample is not present or is below the detection limits of the test.

POSITIVE: Two distinct colored lines appear. One colored line should be in the control region (C) and another colored line should be in the Test region (T). Positive result indicates detection of SARS-CoV-2 antigens in the sample.

Note: The intensity of the color of the test line (T) will vary based on the amount of COVID-19 antigen present in the sample. So any shade of color should be considered positive.

INVALID: If no line appears in the control line region the test should be considered invalid and must be repeated. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

Internal procedural controls are included in the test. A colored line appearing in the control region (C) is an internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

- The COVID-19 Antigen Rapid Test is for professional and in vitro diagnostic use only.
- The COVID-19 Antigen Rapid Test is a screening test for the acute phase of the disease. Samples collected before or after this step may contain antigenic titers below the minimum detection limit of the test and give negative results.
- For optimal test performance, **proper sample collection is critical. Collect the sample correctly and accurately. An inaccurate sample collection can give incorrect (especially negative) results.**
- The accuracy of the test depends on the quality of the sample. False negatives may result from improper sample collection or storage.
- **The sample collection must be carried out by expert and specially trained personnel.**
- Do not add VTM specimen directly to the specimen well of the test cassette.
- This test should be used for detection of SARS-CoV-2 Antigens in human nasopharyngeal specimens as an aid in the diagnosis of patients with suspected SARS-CoV-2 infection in conjunction with clinical presentation. Neither the quantitative value nor the rate of increase in the concentration of SARS-CoV-2 antigens can be determined by this qualitative test.
- If the test result is negative and clinical symptoms persist, it is recommended to re-sample the patient and test again or test with a molecular diagnostic device to rule out infection in these individuals.
- The test will show negative results if the titer of the antigens in the sample is lower than the minimum detection limit of the test.
- In some cases, a positive result may be due to infection with strains of coronaviruses other than SARS-CoV-2.

- **If the result is positive, the test must be confirmed with an alternative method (RT-PCR). Comply with the provisions of the authorities.**
- Excess blood or mucin on the swab specimen may interfere with test performance and may yield a false positive result.
- The COVID-19 Antigen Rapid Test will only indicate the presence of SARS-CoV-2 Antigens in the specimen and should not be used as the sole criteria for the diagnosis of SARS-CoV-2 infections.
- As with all diagnostic tests, all results must be considered with other clinical information available to the physician.

PERFORMANCE

SENSITIVITY, SPECIFICITY E ACCURACY

The COVID-19 Antigen Rapid Test has been compared with the reference RT-PCR. Results showed high sensitivity and specificity.

Method	Results	RT-PCR		Total Results
		Positive	Negative	
Test Rapido COVID-19 Ag	Positive	80	1	81
	Negative	3	120	123
Total Results		83	121	204

Relative Sensitivity: 96,4%

Relative Specificity: 99,2%

Accuracy: 98,0%

SPECIFICITY

The COVID-19 Antigen Rapid Test was tested with the following viral strains. No discernible line at either of the test-line regions was observed at these concentrations:

Description	Concentration
Adenovirus type 3	3.16 x 10 ⁴ TCID50/ml
Adenovirus type 7	1.58 x 10 ⁵ TCID50/ml
Human coronavirus OC43	2.45 x 10 ⁶ LD50/ml
Influenza A H1N1	3.16 x 10 ⁵ TCID50/ml
Influenza A H3N2	1 x 10 ⁵ TCID50/ml
Influenza B	3.16 x 10 ⁶ TCID50/ml
Human Rhinovirus 2	2.81 x 10 ⁴ TCID50/ml
Human Rhinovirus 14	1.58 x 10 ⁶ TCID50/ml
Human Rhinovirus 16	8.89 x 10 ⁶ TCID50/ml
Measles	1.58 x 10 ⁴ TCID50/ml
Mumps	1.58 x 10 ⁴ TCID50/ml
Parainfluenza virus 2	1.58 x 10 ⁷ TCID50/ml
Parainfluenza virus 3	1.58 x 10 ⁸ TCID50/ml
Respiratory syncytial virus	8.89 x 10 ⁴ TCID50/ml

PRECISION

Within-run and Between-run precision has been determined by using three specimens of COVID-19 standard control. Three different lots of COVID-19 Antigen Rapid Test have been tested using negative SARS-COV-2 Antigen weak and SARS-COV-2 Antigen Strong. Ten replicates of each level were tested each day for 3 consecutive days. The specimens were correctly identified >99% of the time.

CROSS-REACTIVITY

The following organisms were tested at 1.0x10⁸ org/ml and all found to be negative when tested with the COVID-19 Antigen Rapid Test:

<i>Arcanobacterium</i>	<i>Pseudomonas aeruginosa</i>
<i>Candida albicans</i>	<i>Staphylococcus aureus subsp. aureus</i>
<i>Corynebacterium</i>	<i>Staphylococcus epidermidis</i>
<i>Escherichia coli</i>	<i>Streptococcus pneumoniae</i>
<i>Moraxella catarrhalis</i>	<i>Streptococcus pyogenes</i>
<i>Neisseria lactamica</i>	<i>Streptococcus salivarius</i>
<i>Nisseria sublava</i>	<i>Streptococcus sp group F</i>

BIBLIOGRAPHY

1. Westgard J.O., Barry P.L., Hunt M.R., Groth T., *A multi-rule Shewhart for quality control in clinical chemistry*, Clinical Chemistry 1981;27:493-501



Comply with EC Directive 98/79/CE

INDEX OF SYMBOLS

4°C	Store between 4-30°C		Use by (YYYY-MM)
	Lot Number		For in vitro diagnostics use only
	Manufacturer		Catalog Number
	Tests per kit		Consult Instruction for Use
	Do not reuse		CE Mark

Rev.3.0.ENG - 2020.10



Delta Chemie Biotechnology sas

Via Nuova Poggioreale, 157 - 80143 - Naples - Italy

TEL/FAX: +39 081 787 1577

E-MAIL: info@deltachemie.it

WEB: www.deltachemie.it

**COMPANY WITH CERTIFIED QUALITY MANAGEMENT SYSTEM
ISO 9001 - ISO 13485**